XP-endo® Retreatment Sequence

1 blister includes/contient/enthält: 1 x (DR1, XP-endo® Shaper 30/.04, XP-endo® Finisher R 30/.00)

- (EN) see www.fkg.ch for additional information
- (FR) informations supplémentaires sur www.fkg.ch
- (DE) weitere Informationen auf www.fkg.ch

XP-endo® Retreatment
Endodontic instruments
Instruments endodontiques
Endodontische Instrumente

FKG REF. 99.X00.10.09A.XX - n°111 - 2019/02

XP-endo® Retreatment

Description - Beschreibung

DR1 

Speed / vitesse / Drehzahl: 800-1000 rpm
Torque / couple / Drehmoment: 1.5 Ncm

 Active tip

Shaper 30/.04 

Speed / vitesse / Drehzahl: 1000-2500 rpm
Torque / couple / Drehmoment: 1 Ncm

Finisher R 30/.00 

Speed / vitesse / Drehzahl: 800-1000 rpm
Torque / couple / Drehmoment: 1 Ncm

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(EN) General information

Indications for use:

These endodontic instruments are intended for use in medical or hospital facilities, by qualified health professionals. This sequence is intended for retreatment.

Contraindications:

Titanium nickel instruments must not be used on individuals with a known allergic sensitivity to these metals.

Adverse reactions: None known.

Warnings and precautions:

- Take 2 x-rays from different angles. A precise appreciation of canal anatomy is very important.
- In multirrooted teeth, begin with the largest canal.
- Single use instruments must not be retreated for subsequent use. That could compromise their integrity.
- Always inspect the instrument(s) before use, and discard (it) them if there are any visible defect(s).
- When an instrument reaches the end of its life, please dispose of it in accordance with the applicable laws and regulations.

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Figure/Abbildung 1

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Figure/Abbildung 2

1. Gently engage DR1 rotating at 800-1000 rpm and 1.5 Ncm torque into the obturation material in order to create a 3-4 mm starting point for the XP-endo Shaper tip (Fig. 1). The active tip of DR1 facilitates the initial penetration.
2. Place one drop of solvent in the coronal space created and wait for at least 1 minute (Fig. 2).

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